

85.12.31

KERT ENVIRONMENTAL SANITATION KERT

concentrate

ASSURE

150 LF

Liquid Germicidal Detergent

Cleans and sanitizes in one application - inhibits growth  
of pathogenic bacteria - low-foaming formulation

GUARANTEE:

n-Alkyl (50% C14, 40% C12, 10% C16) dimethyl benzyl  
ammonium chloride.....15.0%

WARNING

COMMERCIAL

READ THE LABEL BEFORE USING

REG. NO. 16604

P.C.P. ACT

Net Contents

1 Litre

DIN NO. 504408

CORROSIVE

POISON

KERT CHEMICAL INDUSTRIES INC.  
TORONTO, ONTARIO M9L 1M8

Low-foaming ASSURE 150 concentrate is designed expressly to clean and sanitize Century Tubs in hospitals and nursing homes.

Initial Application: Follow the daily procedure but increase the cycle time (Step 3) to between 15 to 20 minutes the first time Tub is cleaned.

#### Daily Procedure

1. Fill Tub with hot water to a level above the upper turbine.
2. After Tub is filled, add entire contents of this bottle (1 litre) to the water in the Tub.
3. Allow Tub to cycle for 10 minutes.
4. Clean sides of Tub, drain and rinse by refilling with fresh water and draining again.

Note: The above daily procedure may have to be repeated more frequently depending on contamination.

NOTICE TO USER - This product is to be used only in accordance with the directions on this label. It is an offence under the Pest Control Products Act to use a control product under unsafe conditions.

PRECAUTION: Do not take internally and avoid contact with eyes.

FIRST AID: If splashed in eyes, flush thoroughly with water for at least 15 minutes and obtain medical attention promptly. If swallowed, drink milk, egg whites, gelatin solution or, if these are not available, large quantities of water. Call a physician immediately.

TOXICOLOGY: Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression and convulsion may be needed.

\*\*\*\*\*

This label transcript service is offered by the Pest Management Regulatory Agency to provide efficient searching for label information. This service and this information do not replace the official hard-copy label. The PMRA does not provide any guarantee or assurance that the information obtained through this service is accurate, current or correct, and is therefore not liable for any loss resulting, directly or indirectly, from reliance upon this service.

+) )