

T1 MOUSE I SOLID

This bait station is resistant to weather and to tampering by children and pets. For use indoors and outdoors.

MOLD RESISTANT • MOISTURE RESISTANT

REGISTRATION NO.: 22135 PEST CONTROL PRODUCTS ACT

DOMESTIC

GUARANTEE:

Diphacinone.....0.005%

* Contains potassium sorbate at 1.00 % and sodium saccharin at 0.01% as preservatives.
Warning, contains the allergens wheat, soy, and peanuts

READ THE LABEL BEFORE USING

(see back panel for additional precautionary statements)

DO NOT ALLOW CHILDREN OR PETS TO HANDLE BAIT STATION

NET CONTENTS:

1 disposable bait station containing 1 x 28 g bait block
up to
16 disposable bait stations containing 16 x 28 g bait blocks

Manufactured by: Bell Laboratories, Inc., 3699 Kinsman Blvd., Madison, WI 53704 USA
Tel: 608-241-0202

PRECAUTIONS: KEEP AWAY FROM CHILDREN, PETS AND LIVESTOCK. Bait station contents may be harmful or fatal if swallowed or absorbed through the skin. Rubber gloves must be worn when collecting and disposing of emptied bait stations and disposing of dead rodents. Avoid contact with eyes, skin or clothing. Wash hands before eating, drinking, chewing gum, using tobacco or using the toilet. Wash skin thoroughly with soap and water after handling. Wash contaminated clothing, separately from other laundry, with soap and water before reuse. **KEEP AWAY FROM FEED AND FOODSTUFFS.**

FIRST AID: If in eyes, hold eye open and rinse slowly and gently with water for 15–20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control centre or doctor for treatment advice. If on skin or clothing, take off contaminated clothing. Rinse skin immediately with plenty of water for 15–20 minutes. Call a poison control centre or doctor for treatment advice. If swallowed, call a poison control centre or doctor immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to do so by a poison control centre or doctor. Do not give anything by mouth to an unconscious person. For all cases of human ingestion, immediately notify a physician or Poison Control Centre. If a pet or livestock poisoning is suspected,

immediately contact a veterinarian. Take container, label or product name and Pest Control Product Registration Number with you when seeking medical attention.

TOXICOLOGICAL INFORMATION: Vitamin K₁ in the form of intramuscular or subcutaneous injections, or by oral ingestion are suggested remedial treatments for anticoagulant poisoning. The severity of the case measured by establishing prolonged prothrombin times (P.T.) will determine appropriate therapy. Monitoring P.T. will indicate the necessity of repeated treatments.

ENVIRONMENTAL HAZARD: This product is toxic to fish and wildlife. Keep out of lakes, streams or ponds.

DIRECTIONS FOR USE

Use bait to control House mouse (*Mus musculus*) indoors and outdoors within 15 metres of buildings.

USE LIMITATIONS: This bait station may be used indoors and outdoors within 15 metres of buildings accessible to children, pets and non-target animals, consistent with all use restrictions and other requirements indicated on this label. **DO NOT** place bait in areas where there is a possibility of contaminating food or surfaces that come in direct contact with food.

SELECTION OF TREATMENT AREAS: Determine areas where mice will most likely find and consume the bait. Generally, these areas are along walls, by gnawed openings, in or beside burrows, in corners and concealed places, between floors and walls, or in locations where rodents or their signs have been observed.

APPLICATION DIRECTIONS

Place the bait station at highest point of activity beside burrow, along a wall or in corner where mice or their signs have been seen. If additional units are to be used, place bait stations at 2.5 to 3.5 metre intervals. Two bait stations may be needed at points of high mouse activity. Maintain an uninterrupted supply of fresh bait for 15 days or until signs of mouse activity cease.

Collect and dispose of unconsumed bait and dead exposed animals according to local requirements. To prevent reinfestation, eliminate food, water, and harborage as much as possible. If reinfestation does occur, repeat treatment. Where a continuous source of infestation is present, establish permanent bait stations and replenish bait as needed.

DISPOSAL: Do not reuse empty container. Wrap unwanted or damaged bait stations and discard in household garbage. **DO NOT** use damaged bait stations.

RODENT DISPOSAL: Dispose of dead rodents in garbage or by burying.

STORAGE: Store in cool, dry place away from other chemicals and food or feed. Store product not in use, in original container, in a secure location inaccessible to children and non-target animals.

NOTICE TO USER: This pest control product is to be used only in accordance with the directions on the label. It is an offence under the *Pest Control Products Act* to use this product in a way that is inconsistent with the directions on the label. The user assumes the risk to persons or property that arises from any such use of this product.

BAIT STATION LABEL

T1 MOUSE I

This Bait Station is Resistant to Weather and to Tampering by Children and Pets. For Use Indoors and Outdoors.

GUARANTEE:

Diphacinone.....0.005%

* Contains potassium sorbate at 1.00 % and sodium saccharin at 0.01% as preservatives.
Warning, contains the allergens wheat, soy, and peanuts

REG. NO. 22135 PEST CONTROL PRODUCTS ACT

Net Contents: 1 x 28g bait block

WARNING



POISON

DO NOT ALLOW CHILDREN OR PETS TO HANDLE BAIT STATION.

Bell Laboratories, Inc., 3699 Kinsman Blvd., Madison, WI 53704 USA
Tel: 608-241-0202

CAUTION : See outer packaging for DIRECTIONS FOR USE, FIRST AID and additional PRECAUTIONARY STATEMENTS.

MODEL NO. MD1

OPTIONAL MARKETING STATEMENTS:

Kills House Mice
Bait station
House Mouse Bait Station
Pre-filled
Pre-baited
Ready to use
Disposable
Pre-baited/prefilled
See thru lid
See thru lid for easy monitoring
With diphacinone

OPTIONAL MARKETING STATEMENTS (CONT'D):

Contains one bait station prefilled with a 28 g bait block
House Mouse killer
Contains 28 g block bait
Pre-filled (pre-baited), ready to use bait station
Feeding station
House Mouse feeding station
Made in USA
Disposable bait station
Can be used indoors or (&) out
Child resistant
Station is child resistant
Station is child and dog resistant
Dog resistant
Child and dog resistant
Child, dog and weather resistant
Weather resistant
Use indoors and out
For use indoors and out
Indoor and outdoor
Indoor and outdoor use
Durable design is child and dog resistant
Kid tested
Dog Tested
Kid and Dog Tested
Kid and Dog Tested Station
Kid and Dog Resistant
Kid and Dog Resistant Station
Kid & Dog Tamper-Resistant
Can be used horizontally or vertically
Clear lid
Clear lid for easy bait monitoring
See-thru lid for easy bait monitoring
No need to touch the bait when placing station

042412