

SpaGuard®

Silk Smart Tabs®

SANITIZER
For Spas & Hot Tubs

Stabilized Chlorinating Tablets
TABLET

Controls bacteria and algae in spa water.
Gentle on equipment. Super convenient to use.
Contains Smart Guard® Chlorine Protector.
Featuring SilkGuard®

GUARANTEE:

Available Chlorine present as Trichloro-s-triazinetriene 85.37%

DOMESTIC

DANGER

POISON

CORROSIVE TO EYES and SKIN

REG. NO. 31250 P.C.P. ACT

READ THE LABEL BEFORE USING

KEEP OUT OF REACH OF CHILDREN

NET CONTENTS: 200 g – 8 kg

KIK HOLDCO COMPANY INC.
33 MacIntosh Blvd.
Concord ON L4K 4L5
1-888-640-7946

DIRECTIONS FOR USE:

Weight per **SpaGuard® Silk Smart Tabs®** tablet: Approximately 20 g

WARNING: Direct **SpaGuard Silk Smart Tabs** contact may cause bleaching. Automatic feeding devices employing **SpaGuard Silk Smart Tabs** must be totally free from traces of other spa sanitizers to avoid the probability of an explosion.

1. Maintain recommended daily levels as determined by testing kit.

Free available chlorine: 3 - 5 ppm

pH: 7.2 - 7.8

Total Alkalinity: 100 - 120 ppm

Calcium Hardness: 200 - 300 ppm

2. Superchlorinate a freshly filled spa by following the manufacturer's instructions of the superchlorinator employed. Maintain water circulation at least 15 minutes before using. Repeat superchlorination procedure after each use period.

3. If using an automatic chlorinator, fill the chamber(s) with **SpaGuard Silk Smart Tabs** and adjust flow rate, according to manufacturer's recommendations, to maintain proper chlorine residual. Adding **SpaGuard Silk Smart Tabs** to the spa's skimmer basket is an alternative chlorine feed method. Daily usage should approximate 10 g of **SpaGuard Silk Smart Tabs** per 1000 L of spa water.

NOTE: Hot weather, heavy spa usage, and use of spa bubblers may require higher usage rates to maintain proper chlorine residuals.

4. Periodic superchlorination of spa water is recommended. Follow label instructions of superchlorinating product being employed.

5. Persistent foam formation, turbid, or malodorous water are indicators of inadequate sanitizer levels. To correct such conditions, superchlorinate spa water by following the manufacturer's instructions of the superchlorinator employed. Maintain circulation at least 30 minutes and wait 4 hours before using. **SpaGuard Silk Smart Tabs** contain a stabilizer (cyanuric acid) which will accumulate in spa water. Keep stabilizer within 30 - 100 ppm range by dilution with fresh water.

6. For proper sanitation, spa must be completely drained periodically. The number of days between COMPLETE SPA DRAINAGE is equal to the volume of spa water in litres, divided by 10 times the maximum number of daily spa users. Refill spa with water and repeat DIRECTIONS FOR USE.

7. Maximum spa water usage temperature is 40°C. Duration in spa water at 40°C should not exceed 15 minutes.

PRECAUTIONS: KEEP OUT OF REACH OF CHILDREN. Dangerous gas formed when mixed with acid. Do not mix with any other chemical. Strong oxidizing agent. May

cause fire or explosion if contaminated. Corrosive to the eye. Severely irritating to the skin. DO NOT get in eyes or on skin. It is recommended to wear rubber gloves when handling this product. Fatal or poisonous if swallowed. Keep container closed and store in a cool, dry location. Do not inhale dust particles.

FIRST AID: **If swallowed**, call a poison control centre or doctor immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to do so by a poison control centre or doctor. Do not give anything by mouth to an unconscious person. **If on skin or clothing**, take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control centre or doctor for treatment advice. **If in eyes**, hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control centre or doctor for treatment advice. Take container, label or product name and Pest Control Product Registration Number with you when seeking medical attention.

TOXICOLOGICAL INFORMATION: Treat symptomatically.

DISPOSAL: Rinse container thoroughly before discarding. Do not reuse empty container. Offer for recycling if available, or dispose in accordance with local and provincial regulations.

For information on disposal of unused, unwanted product, contact the manufacturer or the provincial regulatory agency. Contact the manufacturer and the provincial regulatory agency in case of a spill, and for clean-up of spills.