

2016-3757
2018-03-09

TCCA Alun Spa TABLETS

Swimming pool algicide and bactericide
DOMESTIC

REG. NO. 32984 P.C.P. ACT

KEEP OUT OF REACH OF CHILDREN
READ THE LABEL BEFORE USING

ACTIVE INGREDIENT:

Available chlorine present as trichloro-s-triazinetrione..... 63%

DANGER



POISON



**CORROSIVE
TO EYES, SKIN IRRITANT**

POTENTIAL SKIN SENSITIZER.

Strong oxidant. Contamination or improper use may cause fire, explosion or the release of toxic gases. Add product only to water. DO NOT use in feeders, floaters or any other dispensers that have been used with any other product.

NET CONTENTS: 500g-30 kg

Sani Marc Inc.
42 rue de l'Artisan
Victoriaville, Quebec
G6P 7E3
1-866-470-7264

TCCA Alun Spa contains slow-dissolving stabilized chlorine, clarifying agent and algaecide. The triple action formula of TCCA Alun Spa provides continuous chlorination and effective, long-lasting protection against bacteria and algae growth in spa water. TCCA Alun Spa is specially formulated to help achieve and maintain crystal clear water. TCCA Alun Spa's triple action formula ensures continuous chlorination while clarifying water and providing long-lasting protection against bacteria and algae growth in spa water. Regular use of TCCA Alun Spa will help keep water crystal clear. TCCA Alun Spa eliminates organic matter in spa water such as sunscreen, shampoo, leaves, etc. TCCA Alun Spa's 3-in-1 stabilized chlorine, clarifier and algaecide formula saves on storage space at home and simplifies spa maintenance. The TCCA Alun Spa formula is ideal for clarifying spa water fed by artesian wells. Regular use of TCCA Alun Spa reduces the need for additional clarifying products to keep spa water clear.

For clear spa water, follow our 3-step spa maintenance program: Step 1: Balancing Step 2: Sanitizing and Step 3: Weekly shock treatment

DIRECTI ONS FOR USE

Weight per puck: 15 to 600 g

WARNING:

Direct puck contact may cause bleaching of spa surfaces. Automatic feeding devices employing these pucks must be totally free from traces of other sanitizers to avoid the possibility of an explosion.

1. Maintain recommended daily levels as determined by testing kit.

Free available chlorine: 3.0 – 5.0 ppm

pH: 7.2 – 7.8

Total alkalinity: 100 – 120 ppm

Calcium hardness: 150 – 200 ppm

2. Superchlorinate freshly filled spa by adding 20 g of TCCA Alun Spa per 1,000 L of spa water. Maintain water circulation at least 1 hour before using. (Always ensure free available chlorine is within the recommended levels before using the spa). Repeat superchlorination procedure after each use period.

3. If using an automatic chlorinator, fill the chamber with TCCA Alun Spa pucks and adjust flow rate according to manufacturer's recommendations to maintain proper chlorine residual.

Adding pucks to the skimmer basket is an alternative chlorine feed method.

4. Chlorinate daily to maintain proper chlorine residual with 10 g of TCCA Alun Spa per 1,000 L of spa water, or as needed.

5. Persistent foam formation, turbid, or malodorous water are indicators of inadequate sanitizer levels. To correct such conditions, superchlorinate spa water by adding 30 g of TCCA Alun Spa per 1,000 L of spa water. Maintain water circulation at least 1 hour and wait 4 hours before using. (Always ensure free available chlorine is at the recommended levels before using the spa).

6. Note for superchlorination: TCCA Alun Spa contains a cyanuric acid stabilizer which can accumulate in spa water. Superchlorinating with TCCA Alun Spa will increase the stabilizer level more rapidly. Keep stabilizer level within the range of 30-100 ppm by diluting with fresh water. Superchlorinating with TCCA Alun Spa will also lower total alkalinity. Ensure total alkalinity remains within the recommended levels.

7. NOTE: Maximum spa water usage temperature is 40°C. Duration in spa water at 40°C should not exceed 15 minutes. For proper sanitation, spa must be completely drained periodically. The number of days between COMPLETE SPA DRAINAGE is equal

to the volume of spa water in litres, divided by 10 times the maximum number of daily spa users. Refill spa with water and repeat DIRECTIONS FOR USE.

PRECAUTIONS

KEEP OUT OF REACH OF CHILDREN. Potential skin sensitizer. Strong oxidizing agent. May cause fire or explosion if contaminated. Dangerous gas formed when mixed with acid. Do not mix with any other chemical. It is recommended to wear rubber gloves when handling this product. DO NOT get in eyes, on skin or on clothing. Harmful or fatal if swallowed. Corrosive to eyes. Severely irritating to the skin. This product is irritating to the respiratory tract. Harmful or fatal if inhaled. Avoid inhaling dust and fumes generated by this product.

STORAGE

Keep container closed and store in a cool, dry location.

DISPOSAL

DO NOT reuse the empty containers. Dispose in household garbage. Unused or partially used products should be disposed at provincially or municipally designated hazardous waste disposal sites.

FIRST AID

IF ON EYES: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control centre or doctor for treatment advice.

IF ON SKIN OR CLOTHING: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15 to 20 minutes. Call a poison control centre or doctor for treatment advice.

IF SWALLOWED: Call a poison control centre or doctor immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to do so by a poison control centre or doctor. Do not give anything by mouth to an unconscious person.

IF INHALED: Move person to fresh air. If the person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably mouth-to-mouth, if possible. Call a poison control centre or doctor for treatment advice.

Take container, label or product name and Registration Number with you when seeking medical attention.

TOXICOLOGICAL INFORMATION

Probable mucosal damage may contraindicate gastric lavage. Treat symptomatically.

EASY 2-IN-1

Sanitizes and clarifies spa water
Eliminates algae and helps prevent future growth
Eliminates algae in spa water.
Quickly and continuously eliminates algae in spa water
Destroys algae in spa water quickly and continuously
Provides crystal clear water by destroying bacteria and algae in spa water
Spa algaecide and bactericide
Helps protect against algae growth
Protects chlorine from the sun to extend its effective lifespan
Prevents green water
Kills bacteria in spa water
Long-lasting algaecide
All-in-one Sanitizer, clarifier and Algaecide
All-in-one sanitizer and clarifier
Dual action sanitizer/clarifier
Double action tablets
Double action tablets for routine spa care
Reduces the amount of chlorine needed by preventing organic growth
Leaves no residue
Safe for all spa surfaces
For use in a skimmer or chlorinator
Safe for use in skimmers
Easy to use
Great alternative to liquid sanitizers
Prevents the appearance of black, yellow or mustard algae
Prevents the appearance of black algae
Prevents the appearance of yellow algae
Prevents the appearance of mustard algae
Slow dissolving tablets
Slow dissolving formula makes maintaining chlorine levels easy
Will not stain spa equipment
Better regulates the release of chlorine
Releases chlorine more steadily over time
A stable pH helps to reduce chlorine consumption
Improved chlorine efficiency
Ideal for well water-fed spas
Compatible with salt-based chlorine generators
Compatible with cartridge filters
Improves the efficiency of the filtration system
Quickly resolves cloudy water in spas
Reduces the turbidity of water
Ideal for use with well water
Ensures sparkling clear water

STEP 2 OF AN EASY 3-STEP MAINTENANCE SYSTEM (DRAWING 2)