

2020-0071
2020-04-17

BLUEPRINT SOAK TABS CHLORINE

20 g STABILIZED CHLORINATING TABLETS
CONTROLS BACTERIA AND ALGAE IN SPA WATERS

KEEP OUT OF REACH OF CHILDREN

READ THE LABEL BEFORE USING

CAUTION
POISON
CORROSIVE

KIK HOLDCO COMPANY INC.
33 MacIntosh Blvd.
Concord, ON L4K 4L5
1-888-797-4444

NET CONTENTS
750 g

DOMESTIC
REG. NO. 32351 P.C.P. Act

GUARANTEE:
Available chlorine present as Trichloro-s-triazinetrione ... 90%

PRECAUTIONS: DO NOT MIX WITH OTHER CHEMICALS

Keep out of reach of children. Strong oxidizing agent. May cause fire or explosion if contaminated. Dangerous gas formed when mixed with acid. Do not mix with any other chemical. Avoid contact with skin and eyes. Harmful if swallowed. Keep container closed and store in a cool, dry location. It is recommended to wear rubber gloves when handling this product.

DISPOSAL: Do not re-use empty container. Wrap and dispose of empty container with household garbage. For information on disposal of unused, unwanted product, contact the manufacturer or the provincial regulatory agency. Contact the manufacturer and the provincial regulatory agency in case of a spill, and for clean-up of spills.

FIRST AID DIRECTIONS: In case of contact with eyes or skin, flush thoroughly with water for 15 minutes. For eye contact, get medical attention. If swallowed, drink one glass of milk or water. Never give anything by mouth to an unconscious person. Call a physician or Poison Control Centre IMMEDIATELY. Take container, label or product name and Pest Control Registration Number with you when seeking medical attention.

TOXICOLOGICAL INFORMATION: Treat symptomatically.

Blue Soak Tabs Chlorine are used with either the chlorine or bromine system.

WARNING: Direct tablet contact may cause bleaching of vinyl or painted pool surfaces. Automatic feeding devices employing these tablets must be totally free from traces of other spa sanitizers to avoid the probability of an explosion.

DIRECTIONS FOR USE: Tablet weight: 20 g **WARNING:** Direct tablet contact may cause bleaching of spa surfaces. Automatic feeding devices employing these tablets must be totally free from traces of other spa sanitizers to avoid the probability of an explosion.

Maintain recommended daily levels as determined by testing kit. Free available chlorine: 3-5 ppm, pH: 7.2-7.8 Total alkalinity: 100-120 ppm Calcium hardness: 150-200 ppm. For chlorination, Blueprint Soak Tabs Chlorine should be added to an automatic chlorinator, floating chlorine feeder or spa skimmer to maintain the proper chlorine residual. Superchlorinate a freshly filled spa by adding 20 g of Blueprint Soak Tabs Chlorine per 1,000 L of spa water. Maintain water circulation at least 1 hour before using.

Repeat superchlorination procedure after each use period. Maintain proper chlorine residual by adding 10 g of Blueprint Soak Tabs Chlorine per 1,000 L of spa water, daily or as needed. Persistent foam formation, turbid, or malodorous water are indicators of inadequate sanitizer levels. To correct such conditions, superchlorinate spa water by adding 30 g per 1,000 L spa water. Maintain circulation at least 1 hour and wait four hours before using. Blueprint Soak Tabs Chlorine contains stabilizer (cyanuric acid) which will accumulate in spa water.

Keep stabilizer within 30-100 ppm range by dilution with fresh water. For proper sanitation, spa must be completely drained periodically. The number of days between COMPLETE SPA DRAINAGE is equal to the volume of spa water in litres, divided by 10 times the maximum number of daily spa users. Refill spa with water and repeat DIRECTIONS FOR USE. Note: Maximum spa water usage temperature is 40°C. Duration in spa water at 40°C should not exceed 15 minutes.