

85.12.31

JOHN'S POOL COOL
LIQUID DISINFECTANT FOR
SWIMMING POOL CHLORINATION

WARNING CORROSIVE

NET 20 L

DOMESTIC

REGISTRATION NO. 19011

PEST CONTROL PRODUCTS ACT

GUARANTEE:

Sodium Hypochlorite (when packaged) 10.8%

IMPORTANT: READ LABEL BEFORE USING

JOHN'S RIDGEWAY GENERAL STORE
235 SOUTH MILL STREET
RIDGEWAY, ONTARIO
L0S 1N0
(416) 894-2661

JOHN'S POOL COOL is a high strength sodium hypochlorite solution which exhibits excellent algaecidal properties, proper use of JOHN'S POOL COOL in conjunction with a regular maintenance programme will give a sanitary pool free from harmful micro organisms.

DIRECTIONS

WATER CHLORINATION: Add JOHN'S POOL COOL daily as determined by a test kit to maintain 0.6-1.0 parts per million of residual available chlorine (PPM Cl₂). Daily charge depends upon air temperature and number of people in pool.

TOTAL DAILY CHARGE	LITRES OF WATER
12.5 mL - 25 mL	1,000

125 mL - 250 mL	10,000
250 mL - 500 mL	20,000
625 mL - 1.25 L	50,000
1.25 L - 2.50 L	100,000
3.25 L - 6.25 L	250,000

Wading Pools: Add 5 mL of JOHN'S POOL COOL for every 100 litres of water daily.

Algaecide: Add 375 mL of JOHN'S POOL COOL for every 10,000 L of water and circulate for at least 8 hours. If algae persists, use a weak solution of muriatic acid and scrub walls with brush. Do not use pool until available chlorine has fallen to 2.0 PPM Cl₂.

PH Control: A PH of 7.2 - 7.6 is acceptable. Add muriatic acid if high. Add soda ash if low.

PRECAUTIONS: KEEP OUT OF REACH OF CHILDREN

Contains Sodium Hypochlorite. Avoid contact with skin or eyes. Dangerous gas formed when mixed with acid. Do not mix with ammonia or toilet bowl cleaners, rust remover, household ammonia or acids, and never pre-mix with other pool chemicals. Do not reuse container. Rinse empty container thoroughly with water before disposing.

FIRST AID: If accidental contact occurs flush thoroughly with copious amounts of water. If swallowed, give one or two glasses of milk or warm water, or egg whites. Follow with Emetic and repeat until vomit is clear. Call a physician immediately and keep the patient warm.

STORAGE: Store away from sunlight and in a cool place. USE AT ONCE: Stability maximum 6 months. Do not store in trunk of car for extended periods of time.

This label transcript service is offered by the Pest Management Regulatory Agency to provide efficient searching for label information. This service and this information do not replace the official hard-copy label. The PMRA does not provide any guarantee or assurance that the information obtained through this service is accurate, current or correct, and is therefore not liable for any loss resulting, directly or indirectly, from reliance upon this service.

+))