

93.10.28

LYSOL(R) BRAND

CONCENTRATED

DISINFECTANT

This bottle makes 31 Litres

CLEANS - DISINFECTS - SANITIZES

DOMESTIC

READ THE LABEL BEFORE USING

Keep out of reach of children

REG. NO. 12000 P.C.P. ACT

GUARANTEE:

O-Phenylphenol.....	2.8%
Xylenols.....	1.5%
O-Benzyl-P-Chlorophenol.....	2.7%

350 mL

Distributed by:
L&F Canada Inc.
Aurora, Ontario L4G 3M7

USED FOR GENERATIONS - USE IT EVERY DAY.

For generations Lysol has been a trusted household name. Use Lysol every day for a thorough cleaning.

- NURSERY - diaper pail, changing table, crib, toys
- BATHROOM - toilet bowl, shower, sink
- GENERAL HOUSEHOLD CLEANING - LAUNDRY SANITIZER

DIRECTIONS: CONCENTRATE - MAKES APPROXIMATELY 31 LITRES OF CLEANING SOLUTION. Add 45mL (3 tablespoons) to each 4 litres of water. LAUNDRY SANITIZER. After filling the washer with water

add 250mL (1 cup) to sanitize your wash. Keep out of reach of children.

FIRST AID TREATMENT: DANGER - EYE IRRITANT. If splashed in eyes, immediately flush with water for at least 15 minutes. If swallowed, drink large quantities of milk or water. DO NOT INDUCE VOMITING CALL PHYSICIAN IMMEDIATELY. If on skin wash thoroughly with soap and water.

PHOSPHATE FREE

(R) Reg. Trade Mark

This label transcript service is offered by the Pest Management Regulatory Agency to provide efficient searching for label information. This service and this information do not replace the official hard-copy label. The PMRA does not provide any guarantee or assurance that the information obtained through this service is accurate, current or correct, and is therefore not liable for any loss resulting, directly or indirectly, from reliance upon this service.

+))