

# PRINCIPAL PANEL

## COAST MOUNTAIN

### Spa Tabs

STABILIZED CHLORINATING TABLETS  
CONTROLS BACTERIA AND ALGAE IN SPA WATERS

#### DOMESTIC GUARANTEE:

Trichloro-s-triazinetrione.....100%  
Available chlorine content.....90%  
REG. NO. 23395.06 P.C.P. ACT

NET CONTENTS: 2 kg

Tablet Weight: 20 g

#### PRECAUTIONARY SYMBOLS

#### CAUTION

TRIANGLE

POISON  
(SKULL)

TRIANGLE

CORROSIVE  
(HAND)

KEEP OUT OF REACH OF CHILDREN.  
READ THE LABEL BEFORE USING.

Coast Mountain Spas, #4-3347-262nd St.,  
Langley, B.C. V4W 3V9

## SECONDARY PANEL

### SPA TABS

#### PRECAUTIONS

- KEEP OUT OF REACH OF CHILDREN.
- Strong oxidizing agent. May cause fire or explosion if contaminated.
- Dangerous gas formed when mixed with acid.
- Do not mix with any other chemical.
- Avoid contact with skin and eyes.
- Harmful if swallowed.
- Keep container closed and store in a cool, dry location.

DISPOSAL: Do not re-use empty container. Wrap and dispose of empty container with household garbage. For information on disposal of unused, unwanted product, contact the manufacturer or the provincial regulatory agency. Contact the manufacturer and the provincial regulatory agency in case of a spill, and for clean-up of spills.

FIRST AID: If splashed in eyes or on skin, flush thoroughly with water for 15 minutes. For eye contact, get medical attention. If swallowed, drink 2 or 3 glasses of milk or water. Call a physician or poison control centre IMMEDIATELY. Do not induce vomiting. Take container, label or product name and pest control product registration number with you when seeking medical attention.

TOXICOLOGICAL INFORMATION: Treat symptomatically.

#### DIRECTIONS FOR USE

WARNING: Direct tablet contact may cause bleaching of spa surfaces. Automatic feeding devices employing these tablets must be totally free from traces of other spa sanitizers to avoid the probability of an explosion. Maintain recommended daily levels as determined by testing kit.

Free available chlorine: 3.0-5.0 ppm    pH 7.2-7.8

Total alkalinity: 100-120 ppm    Calcium hardness: 150-200 ppm

For chlorination, SPA-TABS should be added to an automatic chlorinator, floating chlorine feeder or spa skimmer to maintain the proper chlorine residual.

Superchlorinate a freshly filled spa by adding 20 g. of SPA-TABS per 1000 L of spa water. Maintain water circulation at least 1 hour before using. Repeat superchlorination procedure after each use period. Maintain proper chlorine residual by adding 10 g of SPA-TABS per 1000 L of spa water, daily or as needed. Persistent foam formation, turbid, or malodorous water are indicators of inadequate sanitizer levels. To correct such conditions, superchlorinate spa water by adding 30 g of SPA-TABS per 1000 L of spa water. Maintain circulation at least 1 hour and wait four hours before using. SPA-TABS contain stabilizer (cyanuric acid) which will accumulate in spa water. Keep stabilizer within 30 - 100 ppm range by dilution with fresh water. For proper sanitation, spa must be completely drained periodically. The number of days between COMPLETE SPA DRAINAGE is equal to the volume of spa water in litres, divided by 10 times the maximum number of daily spa users. Refill spa with water and repeat DIRECTIONS FOR USE.

Maximum spa water usage temperature is 40 °C. Duration in spa water at 40 °C should not exceed 15 minutes.

\*\*\*\*\*

This label transcript service is offered by the Pest Management Regulatory Agency to provide efficient searching for label information. This service and this information do not replace the official hard-copy label. The PMRA does not provide any guarantee or assurance that the information obtained through this service is accurate, current or correct, and is therefore not liable for any loss resulting, directly or indirectly, from reliance upon this service.