

90.05.29

CLECL - OR
CHLORINATING LIQUID
CONTROLS BACTERIA AND ALGAE
IN SWIMMING POOL WATERS

DOMESTIC

REGISTRATION NO. 21515 PEST CONTROL PRODUCTS ACT

GUARANTEE: (when packaged)

Sodium Hypochlorite.....10.8%
Available Chlorine Content.....10.3%

WARNING

CORROSIVE

IMPORTANT: READ THE LABEL BEFORE USING

Net Contents

CANADIAN LIVESTOCK EQUIPMENT COMPANY LTD.
14 Church St.,
Woodslee, Ont. N0R 1V0
PHONE: (519) 975-2513

DIRECTIONS FOR USE

1. Maintain recommended daily levels as determined by testing kit.

Free available chlorine.....1.0 - 3.0 ppm
pH.....7.2 - 7.8
Total Alkalinity.....100 - 120 ppm
Calcium Hardness.....200 - 300 ppm

2. Superchlorinate pool water initially with 1 L CLECL-OR per 10,000 L or pool water.
3. For best results superchlorinate at dusk.
4. Chlorinate daily to maintain proper chlorine residual with 250 mL of CLECL-OR per 10,000 litres of pool water, or as needed. NOTE: Hot weather, heavy pool usage, and rain may require higher usage rates to maintain proper chlorine residuals.

5. Superchlorinate weekly with 1 L per 10,000 litres when average daily temperature is 25°C - 32 °C. Otherwise, superchlorinate every other week.
6. For visible algae superchlorinate using 2 L per 10,000 L of pool water. Heavy algae infestation may require dosages up to 3L per 10,000 L. Wait at least 24 hours and vacuum pool when dead algae settles.
7. For wading pools, chlorinate daily to maintain chlorine residual of 3-5 ppm by adding 45.0 mL per 1000 L of pool water, or as needed.

NOTE: For outdoor pools, chlorine residuals can be protected from destruction by sun's rays by addition of cyanuric acid (stabilizer).

PRECAUTIONS

- Keep out of reach of children.
- Dangerous gas formed when mixed with acid.
- Do not mix with any other chemical.
- Avoid contact with skin and eyes.
- Harmful if swallowed.
- Keep container closed and store in a cool location.
- Rinse container thoroughly before discarding.
- To dispose of chemical, flush into sewer with large quantities of water.

FIRST AID: If splashed in eyes or on skin, flush thoroughly with water for 15 minutes. For eye contact, get medical attention. If swallowed, drink one glass of milk or water. Call a Physician IMMEDIATELY. Keep patient warm.

STORAGE: Use at once - stability maximum 6 months. Do not store in trunk of car for extended periods of time.

This label transcript service is offered by the Pest Management Regulatory Agency to provide efficient searching for label information. This service and this information do not replace the official hard-copy label. The PMRA does not provide any guarantee or assurance that the information obtained through this service is accurate, current or correct, and is therefore not liable for any loss resulting, directly or indirectly, from reliance upon this service.

+))