

85.12.31

TOP CLEAN
HEAVY-DUTY INSTITUTIONAL STRENGTH
CLEANER - DISINFECTANT - DEODORIZER
FUNGICIDE - VIRUCIDE
AN INLAND ALCARE PRODUCT

GUARANTEE:

Octyl Decyl Dimethyl Ammonium Chloride.....	0.950
Dioctyl Dimethyl Ammonium Chloride.....	0.475
Didecyl Dimethyl Ammonium Chloride.....	0.475

REGISTRATION NO. 16146

PEST CONTROL PRODUCTS ACT

FOR INDUSTRIAL USE ONLY

READ LABEL BEFORE USING

4.55 litre 22.75 litre

MANUFACTURED BY
INLAND ALCARE
10218 - 111 STREET
EDMONTON, ALBERTA T5J 1K9

INLAND/ALCARE TOP CLEAN is a Heavy-Duty Industrial Strength Cleaner formulated to provide maximum cleaning efficiency as it disinfects. It is especially designed for schools, institutions, and industry. A single application removes dirt and grime, disinfects, and is virucidal. When used as directed, the hard surface disinfecting action of TOP CLEAN will reduce the hazard of crossinfection.

When used as directed, TOP CLEAN will deodorize toilet areas, behind and under sinks and counters, garbage cans and garbage storage areas, and places where bacteria growth can cause mal

odors.

PRECAUTIONS:

KEEP OUT OF REACH OF CHILDREN. Harmful if swallowed. May cause eye damage or skin irritation. Do not get concentrate in eyes or on skin. Avoid contamination of food.

FIRST AID: In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes, GET MEDICAL ATTENTION. If swallowed, drink milk, egg whites, gelatin solution, or if these are not available, drink large quantities of water. GET MEDICAL ATTENTION.

NOTICE TO USER:

This control product is to be used only in accordance with the directions on this label. It is an offence under the Pest Control Products Act to use a control product under unsafe conditions.

USE DIRECTIONS

Add 1 part to 40 parts water. Apply TOP CLEAN to walls, floors, and other hard surfaces such as tables, chairs, and bed frames with a cloth or mop.

At 1 part Top Clean to 40 parts water use-level, TOP CLEAN is effective against *Staphylococcus aureus*, *Salmonella choleraesuis*, and *Escherichia coli*. Germicidal performance against these organisms has been confirmed by the AOAC Use-Dilution Test.

NOTE: This product is not effective at the above recommended use-level of four ounces per gallon against *Pseudomonas aeruginosa*.

* At 1 part Top Clean to 40 parts water use-level, TOP CLEAN is virucidal against Herpes Simplex (a member of the virus family that causes infectious mononucleosis), Vaccinia (representative of the pox viruses), and Influenza A2 (the Hong Kong flu Virus), on inanimate environmental surfaces.

At 1 part Top Clean to 22 parts water use-level, TOP CLEAN is effective against *Trichophyton interdigitale*. Fungicidal performance was determined by using the AOAC Fungicidal Test.

Rinse empty container with water before discarding it.

PART OF THE TOP CLEANING SYSTEM

This label transcript service is offered by the Pest Management Regulatory Agency to provide efficient searching for label information. This service and this information do not replace the official hard-copy label. The PMRA does not provide any guarantee or assurance that the information obtained through this service is accurate, current or correct, and is therefore not liable for any loss resulting, directly or indirectly, from reliance upon this service.

+))