

95/2000

PISCI-KILL

CONTROLS BLACK ALGAE

Liquid Algaecide Concentrate

DOMESTIC  
NON FOAMING

READ THE LABEL BEFORE USING

REG. NO.18718.30  
P.C.P. ACT

GUARANTEE:

Poly (oxyethylene (dimethylimino)  
ethylene (dimethylimino) ethylene dichloride)....30%

1 L

MANUFACTURED FOR  
LE MAITRE PISCINIER  
2285, FRANCIS-HUGUES,  
LAVAL (QUEBEC) H7S 1N5  
BY/PAR SANI-MARC INC.

DESCRIPTION:

- Effective on black, green, mustard and other algae.
- Helps flocculate dead algae and sedimentation to pool bottom

- for easy cleaning.
- Effective at any pH level.
- Helps reduce chlorine demand.

DIRECTIONS:

1. Chlorinate pool at least 12 hours before adding PISCI-KILL.
2. For best distribution dilute PISCI-KILL in pail of water and pour around pool.

INITIAL DOSAGE:

1. New/clean pool: 150 mL per 10 000 litres pool water.
2. Algae heavy concentration: 200 mL per 10 000 litres pool water.

WEEKLY MAINTENANCE DOSAGE:

1. Add 50 mL per 10 000 litres pool water.
2. Vacuum to remove dead algae.

CAUTION:

- Keep out of children's reach
- Avoid contact with eyes, skin and clothing
- Avoid contamination of food
- Harmful if swallowed
- Do not pre-mix with other chemicals

FIRST AID: If splashed in eyes or skin flush thoroughly with water. If eye irritation persists seek medical attention. If swallowed, drink promptly a large quantity of milk, egg whites, gelatin solution, or if these are not available drink large quantities of water. Call a physician IMMEDIATELY.

\*\*\*\*\*

This label transcript service is offered by the Pest Management Regulatory Agency to provide efficient searching for label information. This service and this information do not replace the official hard-copy label. The PMRA does not provide any guarantee or assurance that the information obtained through this service is accurate, current or correct, and is therefore not liable for any loss resulting, directly or indirectly, from reliance upon this service.

+))