

2000/2005

SPA TABS

Tablet Weight: 20 gm

STABILIZED CHLORINATING TABLETS

CONTROLS BACTERIA AND ALGAE
IN SPA WATERS

DOMESTIC

REG. NO. 23395.02 P.C.P. ACT

GUARANTEE:

Trichloro-s-triazinetriene	100%
Available Chlorine Content	90%

PEST CONTROL PRODUCTS ACT

4kg

CAUTION

POISON CORROSIVE

READ THE LABEL BEFORE USING

AQUA-TEK
35 B Leeder Ave.,
Coquitlam, B.C.
V3K 3V5

PRECAUTIONS

- Keep out of reach of children.
- Harmful if swallowed.
- Do not mix with any other chemical.
- Avoid contact with skin and eyes.
- Rinse container thoroughly before discarding. Do not reuse.

DISPOSAL

For information on the disposal of unused, unwanted product and on the cleanup of spills, contact the provincial regulatory agency or the manufacturer.

FIRST AID

If splashed in eyes or on skin, flush thoroughly with water for 15 minutes. For eye contact, get medical attention. If swallowed, drink 2 -3 glasses of milk or water. Call a physician or poison control centre IMMEDIATELY. Do not induce vomiting.

DIRECTIONS FOR USE: WARNING: Direct tablet contact may cause bleaching of spa surfaces. Automatic feeding devices employing these tablets must be totally free from traces of other spa sanitizers to avoid the probability of an explosion.

Maintain recommended daily levels as determined by testing kit.

Free available chlorine: 3.0 - 5.0 ppm

pH: 7.2 - 7.8

Total alkalinity: 100 - 120 ppm

Calcium hardness: 150 - 200 ppm

For chlorination, SPA-TABS should be added to an automatic chlorinator, floating chlorine feeder or spa skimmer to maintain the proper chlorine residual. Superchlorinate a freshly filled spa by adding 20 gm of SPA-TABS per 1000 L of spa water. Maintain water circulation at least 1 hour before using. Repeat superchlorination procedure after each use period. Maintain proper chlorine residual by adding 10 g of SPA-TABS per 1000 L of spa water, daily or as needed. Persistent foam formation, turbid, or malodorous water are indicators of inadequate sanitizer levels. To correct such conditions, superchlorinate spa water by adding 30 g of SPA-TABS per 1000 L of spa water. Maintain circulation at least 1 hour and wait four hours before using. SPA-TABS contain stabilizer (cyanuric acid) which will accumulate in spa water. Keep stabilizer within 30 - 100 ppm

range by dilution with fresh water.
For proper sanitation, spa must be completely drained periodically. The number of days between COMPLETE SPA DRAINAGE is equal to the volume of spa water in litres, divided by 10 times the maximum number of daily spa users. Refill spa with water and repeat DIRECTIONS FOR USE. Maximum spa water usage temperature is 40EC. Duration in spa water at 40EC should not exceed 15 minutes.

50553

Aqua-Tek

This label transcript service is offered by the Pest Management Regulatory Agency to provide efficient searching for label information. This service and this information do not replace the official hard-copy label. The PMRA does not provide any guarantee or assurance that the information obtained through this service is accurate, current or correct, and is therefore not liable for any loss resulting, directly or indirectly, from reliance upon this service.

+))