

85.12.31

Duro-Tec  
WOOD PRESERVATIVE  
LIQUID-GREEN

545-277

NET CONTENTS

4,00 L

DOMESTIC

GUARANTEE: 2% COPPER (NAPHTHENATE)

REG. NO. 10321 P.C.P. ACT.

DANGER POISON

Read the label before using

Produced by: LAURENTIDE CHEMICALS INC.  
4650, 12th Avenue,  
Shawinigan-Sud, P.Q.

#### DIRECTIONS FOR USE

The green wood preservative 545-277 is used in a large number of applications to prevent rot and attack of wood by fungi.

545-277 contains a minimum of 2% copper as copper naphthenate.

It is recommended on all wooden surfaces and particularly on dry lumber, beams, sills, docks, boat houses, boats, greenhouses, etc.

Do not mix this product with clear wood preservative.

This product dries in 2 to 3 hours and covers 14,0 to 20,0 litres

per 100 m2.

The green wood preservative 545-277 may be applied by brush, by spray or by immersing the wood for a period of 2 minutes for each cm of penetrations required.

Wash equipment with Chem-Sol 545-153.

PRECAUTIONS: AVOID INGESTION, INHALATION AND SKIN CONTACT. KEEP OUT OF REACH OF CHILDREN.

FIRST AID INSTRUCTIONS: IN CASE OF CONTACT WITH SKIN, WASH WITH RUBBING ALCOHOL. IN CASE OF INGESTION, DO NOT INDUCE VOMITING. CALL A PHYSICIAN IMMEDIATELY OR THE NEAREST POISON CONTROL CENTRE.

TOXICOLOGICAL INFORMATION: TREAT AS IN THE CASE OF AN INGESTION OF KEROSENE.

KEEP AWAY FROM OPEN FLAME OR SPARKS.

\*\*\*\*\*

This label transcript service is offered by the Pest Management Regulatory Agency to provide efficient searching for label information. This service and this information do not replace the official hard-copy label. The PMRA does not provide any guarantee or assurance that the information obtained through this service is accurate, current or correct, and is therefore not liable for any loss resulting, directly or indirectly, from reliance upon this service.

+))