

95/2000

KNOCK-OUT-12

INDUSTRIAL SODIUM HYPOCHLORITE

10 L 20 L

GUARANTEE:

SODIUM HYPOCHLORITE 10.8% WHEN PACKED
(10.3% AVAILABLE CHLORINE CONTENT)

REGISTRATION NUMBER 20573

PEST CONTROL PRODUCTS ACT

READ THE LABEL BEFORE USING

WARNING

CORROSIVE: DANGEROUS GAS FORMED WHEN MIXED
WITH ACID.

MANOTICK POOL & SPA
40 ANN ST.
MANOTICK, ONTARIO
K0A 2N0

HYPOCHLORITE SOLUTION
P.I.N. 1791

USES: For general sanitation, disinfecting, bleaching: use 30 mL of THIS PRODUCT per 5L of water. For specific uses...In water treatment, food and beverage processing, sewage and waste treatment...ask your distributor for technical assistance. FOR SANITARY SWIMMING & WADING POOL WATER: Add the correct amount of THIS PRODUCT to the water daily. THIS PRODUCT is recommended for use in all types of pools: vinyl lined, concrete and fiber glass.

DIRECTIONS FOR USE

REGULAR TREATMENT: Maintain recommended daily levels as determined by testing kit:

Free available chlorine:.....1.0 - 3.0 ppm

pH:.....7.2 - 7.8

Total Alkalinity:.....100 - 120 ppm

Calcium Hardness:.....200 - 300 ppm

1. Before adding THIS PRODUCT test the pH and adjust to within the 7.2 to 7.8 range. Add acid if pH is high. Add soda ash if pH is low.
2. Chlorinate daily to maintain proper chlorine residual with 250 mL of THIS PRODUCT per 10,000 litres of pool water, or as needed. NOTE: Hot weather, heavy pool usage, and rain may require high dosage rates to maintain proper chlorine residuals.

OPENING NEW POOL: Superchlorinate (SHOCK) pool water by adding to the pool 1 L of THIS PRODUCT per 10,000 L of pool water. For best results, superchlorinate (SHOCK) at dusk.

OPENING/CLOSING POOL: For opening/closing pool, superchlorinate (SHOCK) pool water by adding 3 L of THIS PRODUCT per 10,000 L of pool water. For best results superchlorinate (SHOCK) at dusk.

SHOCK TREATMENT: Superchlorinate (SHOCK) using 1 L of THIS PRODUCT per 10,000 litres of pool water when average temperature is 25°C - 32°C. Otherwise, superchlorinate (SHOCK) every other week.

ALGAE CONTROL: For visible algae, superchlorinate (SHOCK) using 3 L of THIS PRODUCT per 10,000 L of pool water. Vacuum pool when dead algae settles.

WADING POOLS: Chlorinate daily to maintain a chlorine residual of 3-5 ppm by adding 45 mL of THIS PRODUCT per 1,000 L of pool water, or as needed.

NOTE: For outdoor pools, chlorine residuals can be protected from destruction by sun's rays by adding a stabilizer (cyanuric acid).

CAUTION: Do not Use THIS PRODUCT on silk, wool, mohair, lycra, unfast colors, leather, aluminum, silverware. Avoid contact with skin and eyes. do not mix with toilet bowl cleaner, rust remover, ammonia or acids.

FIRST AID TREATMENT: contains Sodium Hypochlorite. If splashed in

eyes or on skin, flush thoroughly with water. If swallowed, give one or two glasses of milk or warm water. Call Physician immediately. Keep patient warm.

NOTICE TO USER: This control product is to be used only in accordance with the direction on this label. It is an offence under the Pest Control Products Act to use a control product under unsafe conditions.

STORE UPRIGHT IN A COOL
DARK PLACE

PRODUCT NOT AFFECTED BY
FREEZING

This label transcript service is offered by the Pest Management Regulatory Agency to provide efficient searching for label information. This service and this information do not replace the official hard-copy label. The PMRA does not provide any guarantee or assurance that the information obtained through this service is accurate, current or correct, and is therefore not liable for any loss resulting, directly or indirectly, from reliance upon this service.

+))