

93.06.02

IN THE SWIM(TM)

STICKS

STABILIZED CHLORINATING TABLETS

CONTROLS BACTERIA AND ALGAE
in Swimming Pool Waters

DOMESTIC

REG. NO. 22939, P.C.P. ACT

GUARANTEE:

Trichloro-s-triazinetriene 100%
Available Chlorine Content 90%

CAUTION

POISON CORROSIVE

NET CONTENTS 11.3 Kg

READ THE LABEL BEFORE USING

In The Swim,
320 Industrial Drive
West Chicago, IL
60185

PRECAUTIONS

- Keep out of reach of children.
- Dangerous gas formed when mixed with acid.

- Do not mix with any other chemical.
- Strong oxidizing agent. May cause fire or explosion if contaminated.
- Avoid contact with skin and eyes.
- Harmful if swallowed.
- Keep container closed and store in a cool, dry location.
- Rinse container thoroughly before discarding. Do not reuse.
- To dispose of chemical, flush into sewer with large quantities of water.

FIRST AID

If splashed in eyes or on skin, flush thoroughly with water for 15 minutes. For eye contact, get medical attention. If swallowed, drink one glass of milk or water. Call a Physician IMMEDIATELY.

DIRECTIONS FOR USE

Weight per Tablet 226.80 g

WARNING: Direct Tablet Contact May Cause Bleaching of Vinyl or Painted Pool Surfaces Automatic Feeding Devices Employing These Tablets Must be Totally Free from Traces of other Pool Sanitizers to Avoid the Probability of an Explosion.

1. Maintain recommended daily levels as determined by testing kit.

Free available chlorine:	1.0 - 3.0 ppm
pH:	7.2 - 7.8
Total Alkalinity:	100 - 120 ppm
Calcium Hardness:	200 - 300 ppm

2. For outdoor pools, chlorine residuals can be protected from destruction by the sun's rays by adding 300 g of stabilizer (cyanuric acid) per 10,000 L of pool water to yield the required initial concentration of 30 ppm.
3. If using an automatic chlorinator, fill the chamber(s) with

tablets and adjust flow rate, according to manufacturer's recommendations, to maintain proper chlorine residual. Adding tablets to the pool's skimmer basket is an alternative chlorine feed method. Daily usage should approximate 25 g of Stabilized Chlorinating Sticks per 10,000 L of pool water.

NOTE: Hot weather, heavy pool usage, and rain may require higher usage rates to maintain proper chlorine residuals.

4. Periodic superchlorination of pool water is recommended. Follow label instructions of superchlorination product being employed.
5. STABILIZED CHLORINATING STICKS contain stabilizer which can accumulate in pool water. Keep stabilizer within the range of 30 - 100 ppm by dilution with fresh water.

This label transcript service is offered by the Pest Management Regulatory Agency to provide efficient searching for label information. This service and this information do not replace the official hard-copy label. The PMRA does not provide any guarantee or assurance that the information obtained through this service is accurate, current or correct, and is therefore not liable for any loss resulting, directly or indirectly, from reliance upon this service.

+))