

---

# *Lectranator Digital Model DIG36*

by *AquaCal AutoPilot, Inc.*

**DOMESTIC  
CHLORINE GENERATING DEVICE**  
*Installation and Operation Manual*

MODEL DIG36 REGISTRATION No. 27845 P.C.P. ACT

**Maximum output of hypochlorous equivalent to 0.58 kg of free available chlorine per day.  
One DIG36 Unit Can Treat a Maximum Of 121,000 Litres Of Swimming Pool Water**  
Controls Bacteria And Algae In Swimming Pool Waters

(For Indoor or Outdoor Use)

**RETAIN MANUAL FOR FUTUR REFERENCE**



## ***PRECAUTION:***

***Do not use this device with bromine products.***

***KEEP OUT OF REACH OF CHILDREN***

## **IMPORTANT**

**READ THE LABEL AND THIS MANUAL BEFORE INSTALLING &  
OPERATING**

**Notice to User:** This pest control product is to be used only in accordance with the directions on the label. It is an offence under the Pest Control Products Act to use this product in a way that is inconsistent with the directions on the label. The user assumes the risk to persons or property that arises from any such use of this product.

