

16-OCT-1995

SPA-TABS

STABILIZED CHLORINATING TABLETS

CONTROLS BACTERIA AND ALGAE IN SPA WATERS

REGISTRATION NO. 23395.04

PEST CONTROL PRODUCTS ACT

2 KG

TABLE WEIGHT 20 GM

READ THE LABEL BEFORE USING

CAUTION

POISON CORROSIVE

DOMESTIC

GUARANTEE:

Trichloro-s-triazinetriene	100%
Available Chlorine Content	90%

KOKO BEACH  
THE HOT TUB PEOPLE  
SALES SERVICE

PRECAUTIONS:

- Keep out of reach of children.
- Strong oxidizing agent. May cause fire or explosion if contaminated.
- Dangerous gas formed when mixed with acid.
- Do not mix with any other chemical.
- Avoid contact with skin and eyes.
- Harmful if swallowed.
- Keep container closed and store in a cool, dry location.
- Rinse container thoroughly before discarding. Do not reuse.
- For information on the disposal of unused, unwanted product and the cleanup of spills, contact your municipality, provincial Ministry of Environment or the regional office of Environment Canada.

FIRST AID:

In case of contact in eyes or on skin, flush thoroughly with water for 15 minutes. For eye contact, get medical attention. If swallowed, drink one glass of milk or water. Call a physician or poison control centre IMMEDIATELY.

DIRECTIONS FOR USE:

WARNING: Direct tablet contact may cause bleaching of spa surfaces. Automatic feeding devices employing these tablets must be totally free from traces of other spa sanitizers to avoid the probability of an explosion.

Maintain recommended daily levels as determined by testing kit.

Free available chlorine: 3.0-5.0ppm	Total alkalinity: 100-120ppm
pH 7.2-7.8	Calcium hardness: 150-200ppm

For chlorination, SPA-TABS should be added to an automatic chlorinator, floating chlorine feeder or spa skimmer to maintain the proper chlorine residual. Superchlorinate a freshly filled spa by adding 20 g of SPA-TABS per 1000 L of spa water. Maintain water circulation at least 1 hour before using. Repeat superchlorination procedure after each use period. Maintain proper chlorine residual by adding 10 g of SPA-TABS per 1000 L of spa water, daily or as needed. Persistent foam formation, turbid, or malodorous water are indicators of inadequate sanitizer levels. To correct such conditions, superchlorinate spa water by adding 30 g of SPA-TABS per 1000 L of spa water. Maintain circulation at least 1 hour and wait four hours before using.

SPA-TABS contain stabilizer (cyanuric acid) which will accumulate in spa water. Keep stabilizer within 30 - 100 ppm range by dilution with fresh water. For proper sanitation, spa must be completely drained periodically. The number of days between COMPLETE SPA DRAINAGE is equal to the volume of spa water in litres, divided by 10 times the maximum number of daily spa user. Refill spa with water and repeat DIRECTIONS FOR USE. Maximum spa water usage temperature is 40°C. Duration in spa water at 40°C should not exceed 15 minutes.

Edmonton, Alberta T6A 4Y1

\*\*\*\*\*

This label transcript service is offered by the Pest Management Regulatory Agency to provide efficient searching for label information. This service and this information do not replace the official hard-copy label. The PMRA does not provide any guarantee or assurance that the information obtained through this service is accurate, current or correct, and is therefore not

liable for any loss resulting, directly or indirectly, from  
reliance upon this service.

+) )