

85.12.31

CONCENTRATE ALGICIDE EFFECTIVE
ON STUBBORN ALGAE GROWTH

LIQUID POLYKILL

DOMESTIC

GUARANTEE:

Poly (oxyethylene (dimethylimino) ethylene
(dimethylimino) ethylene dichloride 40%

REGISTRATION NO. 17273
PEST CONTROL PRODUCTS ACT

READ THE LABEL BEFORE USING

1 LITRE

Pool Owners Distributing
Scarborough, Ontario
M1P 2L8

Non-metallic, non-foaming Polykill controls and kills green algae. It's also especially effective for algae which are chlorine resistant and usually cling to pool walls and bottom, such as black and mustard algae.

PRECAUTION: KEEP OUT OF REACH OF CHILDREN. Avoid contact with skin and eyes. Do not take internally. Avoid contamination of food. Do not pre-mix with other chemicals.

FIRST AID TREATMENT: If splashed in eyes or on skin, rinse thoroughly with water. If eye irritation persists, seek medical attention. If swallowed, call physician immediately.

DIRECTIONS: Superchlorinate pool at least 12 hours before applying Polykill. Preventative maintenance: For initial dose, add 75-115 mL per 10,000 litres of pool water. For maintenance dose, add 25-45 ml per 10,000 litres of pool water. Visible algal growth: Add initial dose of 115-190 ml per 10,000 litres of pool water and remove settled algal debris by vacuuming and filtration.

This label transcript service is offered by the Pest Management Regulatory Agency to provide efficient searching for label information. This service and this information do not replace the official hard-copy label. The PMRA does not provide any guarantee or assurance that the information obtained through this service is accurate, current or correct, and is therefore not liable for any loss resulting, directly or indirectly, from reliance upon this service.

+))