

91.10.28

AQUAPOOL
CHLORINATING TABLETS
CONTROLS BACTERIA AND ALGAE IN SWIMMING POOL WATERS

15.9KG

GUARANTEE: Calcium Hypochlorite 70%
Available Chlorine Content 70%

WARNING

READ THE LABEL BEFORE USING

POISON

CORROSIVE

DOMESTIC REG. NO. 13753 P.C.P. ACT

Canadian Miraclean Products Ltd.
12431 Horseshoe Way
Richmond, British Columbia V7A 4X6

PRECAUTIONS: Keep out of reach of children.

Dangerous gas formed when mixed with acid.

Do not mix with any other chemical.
Strong oxidizing agent. May cause fire or explosion if
contaminated.

Avoid contact with skin and eyes- harmful if swallowed.

Keep container closed and store in a cool and dry location away
from heat or sunlight.

Rinse container thoroughly before discarding. Do NOT reuse
container. To dispose of chemical, flush into sewer with large
quantities of water.

First Aid: If splashed in eyes or on skin, flush thoroughly
with water for 15 minutes. For eye contact, get medical
attention. If swallowed, drink two or three glasses of milk or
water. Call a physician immediately. Keep patient warm.

DIRECTIONS FOR USE: Weight per tablet: 20 grams.

WARNING: Direct tablet contact may cause bleaching of vinyl or painted pool surfaces.

Automatic feeding devices employing these tablets must be totally free from traces of other pool sanitizers to avoid probability of explosion.

1. Maintain proper daily levels as determined by testing kit.

Free available chlorine:	1.0 - 3.0 ppm
pH:	7.2 - 7.8 ppm
Total Alkalinity:	80 - 120 ppm
Calcium Hardness:	200 - 300 ppm

2. To minimize bleaching of vinyl liners, predissolve required amount of Aqua in a plastic pail with concentration not exceeding 1kg per 10L. Always add chemical to water, never water to chemical. Dissolution time is approximately 4 hours. Add resulting slurry directly to pool water.
3. Initially, superchlorinate pool water by adding 150g of the predissolved Aqua per 10,000L of pool water. For best results superchlorinate at dusk.
4. Maintain proper chlorine residual by broadcasting daily over the pool's surface. 30 g of the predissolved Aqua Tablets per 10,000L of pool water, or as needed. For automatic feeding devices, adjust flow rate, according to manufacture's recommendations, to maintain proper chlorine residual.

Note: Hot weather, heavy pool usage, and rain may require higher usage rates to maintain proper chlorine residuals.

5. Superchlorinate weekly with 150g of the predissolved tablets per 10,000L when average daily temperature is 25'C - 32'C, otherwise, superchlorinate every other week.
6. For visible algae superchlorinate using 350g of the predissolved tablets per 10,000L of pool water. Heavy algae infestation may require dosages up to 500g per 10,000L. Vacuum pool when dead algae settles.

Note: For outdoor pools, chlorine residuals can be protected from destruction by sun's rays by addition of stabilizer (cyanuric acid).

A PRODUCT OF
CANADIAN MIRACLEAN PRODUCTS LTD.
12431 HORSESHOE WAY,
RICHMOND BRITISH COLUMBIA
V7A 4X6

This label transcript service is offered by the Pest Management Regulatory Agency to provide efficient searching for label information. This service and this information do not replace the official hard-copy label. The PMRA does not provide any guarantee or assurance that the information obtained through this service is accurate, current or correct, and is therefore not liable for any loss resulting, directly or indirectly, from reliance upon this service.

+))