

91.06.11

BINBROOK CHLORINATING LIQUID

Controls Bacteria and Algae in Swimming Pool Waters

DOMESTIC

REGISTRATION NO. 21140  
PEST CONTROL PRODUCTS ACT

Net Contents:

IMPORTANT: READ THE LABEL BEFORE USING!

WARNING CORROSIVE

Made in Canada

BINBROOK SUPPLIES  
3194 Hwy. 56, Binbrook, Ontario  
L0R 1C0

Telephone: 692-3535

BINBROOK CHLORINATING LIQUID: is a high strength sodium hypochlorite solution which exhibits excellent algaecidal properties. Proper use of BINBROOK CHLORINATING LIQUID in conjunction with a regular maintenance program will give a sanitary pool, free of harmful micro-organisms.

DIRECTIONS

1. MAINTAIN PROPER DAILY LEVELS AS DETERMINED BY TESTING KIT:

FREE AVAILABLE CHLORINE.....	1.0 - 3.0 ppm
pH .....	7.2 - 7.8
TOTAL ALKALINITY .....	100 - 200 ppm
CALCIUM HARDNESS .....	200 - 300 ppm

2. SUPERCHLORINATE POOL WATER INITIALLY WITH 1 LITRE OF BINBROOK CHLORINATING LIQUID PER 10,000 LITRES OF POOL WATER.
3. FOR BEST RESULTS SUPERCHLORINATE AT DUSK.

4. CHLORINATE DAILY TO MAINTAIN PROPER CHLORINE RESIDUAL WITH 250 MILLILITRES OF BINBROOK CHLORINATING LIQUID PER 10,000 LITRES OF POOL WATER, OR AS NEEDED. NOTE: HOT WEATHER, HEAVY POOL USAGE AND RAIN, MAY REQUIRE HIGHER DOSAGE RATES TO MAINTAIN PROPER CHLORINE RESIDUALS.
5. SUPERCHLORINATE WEEKLY WITH 1 LITRE PER 10,000 LITRES WHEN AVERAGE DAILY TEMPERATURE IS 25'C - 32'C. OTHERWISE, SUPERCHLORINATE EVERY OTHER WEEK.
6. FOR VISIBLE ALGAE, SUPERCHLORINATE WITH 2 LITRES PER 10,000 LITRES OF POOL WATER. HEAVY ALGAE INFESTATION MAY REQUIRE DOSAGES UP TO 3 L PER 10,000 L. WAIT AT LEAST 24 HOURS AND VACUUM POOL WHEN DEAD ALGAE SETTLES.
7. FOR WADING POOLS, CHLORINATE DAILY TO MAINTAIN CHLORINE RESIDUAL AT 3-5 PPM BY ADDING 45 ML PER 1,000 L OF POOL WATER, OR AS NEEDED.

NOTE:

FOR OUTDOOR POOLS, CHLORINE RESIDUALS CAN BE PROTECTING FROM DESTRUCTION BY SUN'S RAYS BY ADDITION OF CYANURIC ACID STABILIZER.

PRECAUTIONS:

- KEEP OUT OF REACH OF CHILDREN.
- DANGEROUS GAS FORMED WHEN MIXED WITH ACID.
- AVOID CONTACT WITH SKIN AND EYES.
- HARMFUL IF SWALLOWED.
- KEEP CONTAINER CLOSED AND STORE IN A COOL LOCATION.
- RINSE CONTAINER THOROUGHLY BEFORE DISCARDING.
- TO DISPOSE OF CHEMICAL, FLUSH INTO SEWER WITH LARGE QUANTITIES OF WATER.

FIRST AID:

IF SPLASHED IN EYES OR ON SKIN, FLUSH THOROUGHLY WITH WATER FOR 15 MINUTES. FOR EYE CONTACT GET MEDICAL ATTENTION. IF SWALLOWED, DRINK ONE OR TWO GLASSES OF MILK OR WARM WATER. CALL A PHYSICIAN IMMEDIATELY. KEEP PATIENT WARM.

STORAGE: USE AT ONCE - STABILITY MAXIMUM 6 MONTHS. DO NOT STORE IN TRUNK OF CAR FOR EXTENDED PERIODS OF TIME.

GUARANTEE:  
(when packaged)

SODIUM HYPOCHLORITE .....	10.8%
AVAILABLE CHLORINE CONTENT .....	10.3%

\*\*\*\*\*

This label transcript service is offered by the Pest Management Regulatory Agency to provide efficient searching for label information. This service and this information do not replace the official hard-copy label. The PMRA does not provide any guarantee or assurance that the information obtained through this service is accurate, current or correct, and is therefore not liable for any loss resulting, directly or indirectly, from reliance upon this service.

+))